Meeting of the ERFP Ex situ Conservation WG Zagreb, 23-26.05. 2012

Introduction to CryoWEB and state of implementation in Europe

Zhivko Duchev

National gene banks

- Choosing breeds for conservation
- Organization of the material collection
- Freezing and storing the samples
- Documentation system

What is CryoWEB

Open source web-software intended to serve as a electronic register of the national collections, uniformly applicable across various species, types of genetic material, storage facilities.

Develped in Institute of Farm Animal Genetics (FLI) in Germany

With CryoWEB the manager can

- Find the distribution of each recorded sample in the storage facilities
- Get all the information about the sample by the label on its physical storage vessel
- Find the procedure to be followed when thawing the sample
- Perform complete or partial inventory of the storage

With CryoWEB the manager can

- Identify all samples from a single donor
- Keep track of material relocation and use
- Virtually browse the storage facilities
- Keep contact data of involved organizations

Other CryoWEB features

- The page can be accessed from anywhere via Internet
- Access can be restricted to Intranet
- Data access control
- Can be translated and used in the local language

Data organization in CryoWEB

- Core "passport" data
- Extended normalized data
- Extended archived data

Core data

- Essential for the gene bank operation
- Uniform across species and breeds
- Minimal data set

Examples: sample id, donor, vessel type, location

Extended normalized data

- Non essential for the gene bank operation
- Not available for all donors
- Useful additional information


Examples: donor's sire and dam, birth date, coordinates of the donor location

Extended archived data

- Intended for storage only
- Not directly processed in the system
- Organized as free content fields, archives, images

Examples: donor's scaled image, archived documents


Animal form




Animal form II

ANIMAL 0 00 G-003340 Sire ID G-000250 Dam ID G-097797 Animal ID* ♦ ? ♦ ? Sheep RPL Species* Breed* 26-03-2003 Birthday Birthyear \$? ♦ ? Sex* male Organization* ZV-NS-Schaf 52.847123 9.12723 Longitude Latitude Photo (?) Browse... Comments File (?) Browse... View


Basic blocks of mandatory data




CryoWEB snapshots



CryoWEB snapshots II



CryoWEB project page



http://cryoweb.tzv.fal.de/

Additional information

CryoWEB – User's Guide and Reference manual

First published 2010


In Series: Applications in Biodiversity Informatics Institute of Farm Animal Genetics (FLI), Mariensee ISBN 978-3-9813280-0-4

CryoWEB demo site http://cryo-devel.tzv.fal.de/

Implementation in Europe

- Within EFABISnet project
- Outside EFABISnet project

EFABISnet project



http://efabisnet.tzv.fal.de/

Installations in EFABISnet

- Austria
- Estonia
- Finland
- Georgia
- Iceland


- Ireland
- Netherlands
- Slovakia
- Slovenia
- Switzerland

Additional installations

- Germany
- Greece
- Croatia

- Cyprus
- Italy
- Poland

Map of installations in Europe



Use cases for CryoWEB

- Countries having operational gene banks
 - Data migration
 - Data verification
- Countries just starting
 - Model for data collection
 - Ready solution for documentation
- Sub-regional genebanks
 - Italian animal genetic resources cryobank network
 - MTT Finland

Thank you for your attention!